

PRISTINE[®]

High Purity Solvents

LC-MS Reagents

No LC-MS TIC signals higher than 100ppb Reserpine

No LC-MS TIC signals higher than:


50ppb Reserpine (ESI + mode)

50ppb 4-Nitrophenol (ESI – mode)

Very low metal concentration (<100ppb) and low particles

Low Particles

Actual Lot Analysis printed on each label



Part # 1934-5
Isopropyl Alcohol, LC-MS Grade
4 Liter, 4/case

 U.S. roots.
Global routes.