Birch® BIOTECH

# PRISTINE®

## High Purity Solvents

## LC-MS Reagents

- No LC-MS TIC signals higher than 100ppb Reserpine
- No LC-MS TIC signals higher than:
  - 50ppb Reserpine (ESI + mode)
  - 50ppb 4-Nitrophenol (ESI – mode)
- Very low metal concentration (<100ppb) and low particles
- Low Particles
- Actual Lot Analysis printed on each label




Part # 1935-5
Methanol, LC-MS Grade
4 Liter, 4/case

U.S. roots.
Global routes.